

Lipo Surg



400 ml/min infiltration

up to 4200 strokes

Power Liposuction


LipoSurg

Vibration liposuction with tumescent infiltration

The LipoSurg system combines the functions of infiltration and motor-assisted vibration liposuction into a single, compact device. The desired values for infiltration and liposuction are set by selecting the function and pressing the ascending or descending arrow keys. Once the values are set, the liposuction or infiltration can be controlled by the VARIO or the ON/OFF pedal.

The extraction of fat requires the connection of this compact system with the Nouvag high vacuum pump "Vacuson 60". The high flow rate of 60 liters per minute is ideal for efficient liposuction.

Integrated high volume infiltration system with balanced pressure for tissue protection and a flow rate of 400 ml/min.



Lipo Cannulae Range



CETTA

PRIMO



CLOVERLEAF MEDICAL

Exceeding Expectation

Cloverleaf is the UK's leading supplier of Liposuction equipment and instrumentation, we are also very proud of the trading partnerships forged in our ever expanding export market.

A family company that believes strongly in the quality of our products and the longevity of our relationships, our aim is to exceed your expectations, and you can trust our advice as we have over 100 years of combined experience in this speciality alone.

CETTA

The outstanding quality of the CETTA disposable product range is immediately apparent; high grade stainless steel, reassuring weight, excellent balance and superb finish. We strongly believe CETTA is not only the best disposable product available worldwide but is also far superior to the majority of re-usable brands. For Cetta products please refer to the product codes in blue.

PRIMO

The reusable Primo range is as versatile as you need it to be, we will create any bespoke size and style upon request. The range shown in this catalogue is a guide to the most popular products we currently provide, many are regular stock items for fast delivery. Durable, high quality stainless steel shafts, reassuringly weighty and balanced Aluminium handles. For Primo products please refer to the products in gold.

FLUID INFUSION

Infusion needles

High grade stainless steel. 6 spaced holes in a spiral around the tip. For fluid infusion using leur lock syringes, handles or giving set tubing.



CETTA RANGE		PRIMO RANGE	
Product Code	Product Code	Diameter (mm)	Length (cm)
CETTA INF/N 215	PRIMO INF/N 215	2	15
CETTA INF/N 223	PRIMO INF/N 223	2	23
CETTA INF/N 230	PRIMO INF/N 230	2	30
CETTA INF/N 315	PRIMO INF/N 315	3	15
CETTA INF/N 323	PRIMO INF/N 323	3	23
CETTA INF/N 330	PRIMO INF/N 330	3	30

Curved



Product Code	Product Code	Diameter (mm)	Length (cm)
CETTA INF/N1.623CRV	PRIMO INF/N1.623CRV	1.6	23

Infusion handle non valve

Comfortable handle to link tubing to needle, lack of valve reduces leakage. Single shaft allows infiltration and suction.

Product Code
PRIMO INF/HAN W/O



Primo handle: Aluminium

Product Code
CETTA INF/HAN W/O




Cetta handle: POM

Infusion handle with Valve

This infusion handle, with on/off button, offers maximum control for infusion during the tumescent technique.


Product Code
PRIMO HUN/HAN



POWERED SUCTION

One piece cannulae

High grade Stainless Steel shaft, for aspiration with powered suction.



Primo handle: Aluminium

CETTA RANGE		PRIMO RANGE	
Product Code	Product Code	Diameter (mm)	Length (cm)
CETTA 210MM	PRIMO 210MM	2	10
CETTA 315MM	PRIMO 315MM	3	15
CETTA 325MM	PRIMO 325MM	3	25
CETTA 430MM	PRIMO 430MM	4	30
CETTA 530MM	PRIMO 530MM	5	30

Product Code	Product Code	Diameter (mm)	Length (cm)
CETTA 325/30holeM	PRIMO 325/30holeM	3	25

Product Code	Product Code	Diameter (mm)	Length (cm)
CETTA 325LM	PRIMO 325LM	3	25
CETTA 430LM	PRIMO 430LM	4	30
CETTA 530LM	PRIMO 530LM	5	30

Product Code	Product Code	Diameter (mm)	Length (cm)
CETTA 210TSM	PRIMO 210TSM	2	10
CETTA 320TSM	PRIMO 320TSM	3	20
CETTA 430TSM	PRIMO 430TSM	4	30

Product Code	Product Code	Diameter (mm)	Length (cm)
CETTA 210SM	PRIMO 210SM	2	10
CETTA 315SM	PRIMO 315SM	3	15

Product Code	Product Code	Diameter (mm)	Length (cm)
CETTA 215OHM	PRIMO 215OHM	2	15
CETTA 320OHM	PRIMO 320OHM	3	20
CETTA 430OHM	PRIMO 430OHM	4	30

Mercedes

3 holes spaced evenly around the tip. The most popular tip style.



30 Hole harvest cannula

Harvest fat in smaller parcels, reduces blocking.



Lantern

Chinese lantern shaped tip, useful for fibrous areas and popular in Gynecomastia procedures.



The Scream

Two holes above a single hole, a cannula lookalike of the famous Munch painting.



Submental

One hole, flat backed for areas closer to the skin.



One hole

A single downward facing aperture.



TOOMEY SYRINGE SUCTION



Toomey Cannulae
High grade Stainless steel shaft,
for aspiration with 60cc Toomey syringes.

Cetta handle: POM

Primo handle: Aluminium

Mercedes tip

3 holes spaced evenly around the tip. The most popular tip style.



CETTA RANGE	PRIMO RANGE		
Product Code	Product Code	Diameter (mm)	Length (cm)
CETTA 210MTS	PRIMO 210MTS	2	10
CETTA 315MTS	PRIMO 315MTS	3	15
CETTA 318MTS	PRIMO 318MTS	3	18
CETTA 325MTS	PRIMO 325MTS	3	25
CETTA 430MTS	PRIMO 430MTS	4	30

30 Hole harvest cannula

Harvest fat in smaller parcels, reduces blocking.



Product Code	Product Code	Diameter (mm)	Length (cm)
CETTA 2.515/30holeTS	PRIMO2.515/30holeTS	2.5	15
CETTA 320/30holeTS	PRIMO320/30holeTS	3	20

Submental

One hole, flat backed for areas closer to the skin



Product Code	Product Code	Diameter (mm)	Length (cm)
CETTA 210SMTS	PRIMO 210SMTS	2	10
CETTA 315SMTS	PRIMO 315SMTS	3	15

Open harvest

Guarded fully open tip, designed specifically for tissue harvest.



SYRINGE LOCKS

All 60cc locks also fit 50cc BD

Square syringe locks

Product Code	Syringe Size(cc)
PRIMOSQUARE10	10cc
PRIMOSQUARE30	30cc
PRIMOSQUARE60	60cc



60cc syringe lock - single arm

Product Code	Product Code	Syringe Size(cc)
CETTA SNAP10	PRIMO SNAP10	10cc
CETTA SNAP20	PRIMO SNAP20	20cc
CETTA SNAP60	PRIMO SNAP60	60cc


60cc syringe lock - double arm

Fixes to syringe base for added security.

Product Code	Product Code	Syringe Size(cc)
CETTA DOUBSNAP60	PRIMO DOUBSNAP60	60cc

60cc syringe lock - 2 step

Product Code	Product Code	Syringe Size(cc)
CETTA2STEP	PRIMO2STEP	60cc



MICRO HARVEST

Tissue harvest using luer lock syringes.
All Cloverleaf luer lock bases have a thumb platform for comfort and orientation.

Open harvest

Guarded fully open tip, designed specifically for tissue harvest.

CETTA RANGE	PRIMO RANGE		
Product Code	Product Code	Diameter (mm)	Length (cm)
CETTA MIC315/HARV	PRIMO MIC315/HARV	3	15
CETTA MIC322/HARV	PRIMO MIC322/HARV	3	22



30 Hole harvest cannula

Harvest fat in micro droplets for fine graft re-injection

Product Code	Product Code	Diameter (mm)	Length (cm)
CETTA 2.515LL/30hole	PRIMO 2.515LL/30hole	2.5	15
CETTA 315LL/30hole	PRIMO 315LL/30hole	3	15

Submental

One hole, flat backed for areas closer to the skin

Product Code	Product Code	Diameter (mm)	Length (cm)
CETTA 322S/LL	PRIMO 322S/LL	3	22



Mercedes tip

3 holes spaced evenly around the tip.
The most popular tip style.

Product Code	Product Code	Diameter (mm)	Length (cm)
CETTA 215MLL	PRIMO 215MLL	2	15
CETTA 315MLL	PRIMO 315MLL	3	15
CETTA 322MLL	PRIMO 322MLL	3	22



CELL SELECT CANNULAE RANGE

Level 1 - Gentle Rasp

Product Code	Product Code	Holes	Diameter (mm)	Length (cm)
CETTA 2515G	PRIMO 2515G	21	2.5	15
CETTA 320G	PRIMO 320G	30	3	20



Level 2 - Medium Rasp

Product Code	Product Code	Holes	Diameter (mm)	Length (cm)
CETTA 2515M	PRIMO 2515M	21	2.5	15
CETTA 320M	PRIMO 320M	30	3	20

Level 3 - Aggressive Rasp

Product Code	Product Code	Holes	Diameter (mm)	Length (cm)
CETTA 2515A	PRIMO 2515A	21	2.5	15
CETTA 320A	PRIMO 320A	30	3	20


The ingenious Cell Select micro harvest cannula range is specifically designed to capture smaller aliquots of healthy fat cells via numerous 1mm micro-ports. The harvest of smaller cell packages results in easier, smoother re-injection through finer injection cannula. The spacing and increased surface area of the apertures per cannula tip results in an effective and constant collection speed, and greatly reduces blocking.



MICRO INFILTRATION CANNULAE

High grade stainless steel for re-injection using Luer lock syringes.

CETTA RANGE		PRIMO RANGE	
Product Code	Product Code	Diameter (G)	Length (cm)
CETTA MIC16/A7	PRIMO MIC16/A7	16G	7
CETTA MIC16/A9	PRIMO MIC16/A9	16G	9
CETTA MIC16/B7	PRIMO MIC16/B7	16G	7
CETTA MIC16/B9	PRIMO MIC16/B9	16G	9
CETTA MIC16/B12	PRIMO MIC 16/B12	16G	12
CETTA MIC16/C7	PRIMO MIC16/C7	16G	7
CETTA MIC16/C9	PRIMO MIC16/C9	16G	9
CETTA MIC18/A9	PRIMO MIC18/A9	18G	9
CETTA MIC18/B9	PRIMO MIC18/B9	18G	9



Curved Cannula

Product Code	Product Code	Diameter (G)	Length (cm)
CETTA CURV/18B7	PRIMO CURV/18B7	18G	7
CETTA CURV/16B12	PRIMO CURV/16B12	16G	12
CETTA CURV/16B15	PRIMO CURV/16B15	16G	15

Ultra fine

Ultra fine facial needle with unblocking Stylet


Product Code	Product Code	Diameter (G)	Length (cm)
CETTA 235FN/LL	PRIMO 235FN/LL	23G (0.7mm)	5
CETTA 217FN/LL	PRIMO 217FN/LL	21G (0.8mm)	7
CETTA 2185FN/LL	PRIMO 2185FN/LL	21G (0.8mm)	8.5
CETTA 15FN/LL	PRIMO 15FN/LL	20G (1.0mm)	5
CETTA 207FN/LL	PRIMO 207FN/LL	20G (1.0mm)	7
CETTA 189FN/LL	PRIMO 189FN/LL	18G (1.2mm)	9



Micro Drop

Curved Obtuse tip for accurate placement of micro droplets.

Product Code	Product Code	Diameter (mm)	Length (cm)
CETTA 1285MD	PRIMO 1285MD	1.2	8.5



PRIMO or CETTA we can produce excellent products quickly so...
Why not design your own **unique** cannula

DISSECTORS

V-Dissector

Undermining fibrous tissue, widens at the tip and ideal for breast surgery.

CETTA RANGE		PRIMO RANGE	
Product Code	Product Code	Diameter (G)	Length (cm)
CETTA VDISS167	PRIMO VDISS167	16G	7
CETTA VDISS169	PRIMO VDISS169	16G	9

W-Dissector

For undermining small areas of fibrous tissue, less aggressive and finer than V Dissector.

Product Code	Product Code	Diameter (G)	Length (cm)
CETTA WDISS	PRIMO WDISS	16G	8.5

HAY Tunnel Dissector

For blunt dissection tunneling when a V-dissector may be considered aggressive i.e. around Implants. Slight upward curve to tip in order to inhibit direct forceful contact, comfortable handle with orientation marking.

Product Code	Product Code	Diameter (mm)	Length (cm)
CETTA HTD	PRIMO HTD	2	12



ADAPTORS

Toomey to luer lock adaptor

Transfer aspirate from 60cc Toomey to any size luer lock syringe.

Product Code	Product Code
CETTA MIC/LL/S	PRIMO MIC/LL/S



Luer lock to luer lock adaptor

Transfer between any size luer lock syringes.

Product Code	Product Code
CETTA MIC/LL/LL	PRIMO MIC/LL/LL



Toomey syringe to luer lock needle adaptor

Product Code	Product Code
CETTA MIC/LL/N	PRIMO MIC/LL/N



ACCESSORIES

Trigger injection gun 10cc/5cc

Trigger mechanism controlled injection from 10cc BD syringe or 5cc BD syringe (adaptor included), eases the physical pressures of re-injection and controls accidental over filling.


Product Code	Capacity (cc)
PRIMO TG105	5, 10



Syringe rack

Autoclavable steel, holds 60 10cc syringes or 12 60cc syringes.

Product Code
PRIMO MIC/10CCR
PRIMO MIC/60CCR



The golden gun 10cc

Enjoy being Scaramanga with your very own golden gun. Thumb controller injection from 10cc BD syringe, eases the physical pressures of re-injection and controls accidental over filling.

Product Code	Capacity (cc)
PRIMO GG10	10cc



NANOFAT AND MICROFAT SYSTEM



Nanofat and Microfat System es un sistema reutilizable especialmente diseñado para las técnicas nanofat y microfat. Consta de tres conectores de trasvase de 1.2 mm., 1.0 mm. y 0.8 mm. para realizar la técnica microfat y reducir el tamaño de los adipocitos. Incorpora también un dispositivo de filtrado con malla para eliminar las fibras de colágeno y obtener un implante nanofat.

LAS DOS TÉCNICAS EN UN SOLO SISTEMA.

REFERENCIA: PRIMO/NMS